

# Micro-scale flow cytometry-based and biochemical analysis of lipid signaling in primary B cell subpopulations

Pierre Antony<sup>1\*</sup>, Kristen Hoek<sup>1\*</sup>, Bhaskarjyoti Sarmah<sup>2</sup> and Wasif N. Khan<sup>1#</sup>

\*These authors contributed equally.

<sup>1</sup>Department of Microbiology and Immunology and <sup>2</sup>Department of Cell and Developmental Biology, Vanderbilt University School of Medicine, Nashville, TN 37232-0146, USA.

#Corresponding Author: Wasif N. Khan, PhD., Department of Microbiology and Immunology, Vanderbilt University School of Medicine, Nashville, Tennessee 37232-0146, USA. Phone: (615) 343 5632; Fax : (615) 343 7392; E-mail: [wasif.khan@vanderbilt.edu](mailto:wasif.khan@vanderbilt.edu)

Submitted: May 29, 2007; Revised: October 25, 2007; Accepted: December 9, 2007

Indexing terms: receptors, antigen, b-cell; immunomagnetic separation

Abbreviations: BCR, B cell antigen receptor; DAG, diacylglycerol; DMEM, Dulbecco's modified Eagle's media; IP<sub>3</sub>, inositol triphosphate; FCM, flow cytometry; MACS, magnetic bead cell sorting; MFI, mean fluorescence intensity; PKC, protein kinase C; RBC, red blood cell

---

## ABSTRACT

B cell subpopulations in the spleen have been extensively characterized phenotypically; however, biochemical properties of these cell populations following B cell antigen receptor engagement have not been fully determined due to technical difficulties and limiting cell numbers. We therefore employed mini-scale protocols to assess lipid signaling, particularly that of diacylglycerol and inositol triphosphate, with as few as  $0.5 \times 10^6$  purified early (T1) and late (T2) transitional B cells. Additionally, utilizing flow cytometric techniques, we determined levels of phosphatidylinositol bisphosphate and calcium mobilization in T1 and T2 cells, as well as mature follicular and marginal zone B cells using less than  $1 \times 10^6$  primary B cells. Thus, these biochemical and flow cytometric methodologies can be used to analyse signal-induced changes in phosphatidylinositol bisphosphate levels, diacylglycerol and inositol triphosphate production and calcium in each B cell population.

---

## INTRODUCTION

The B cell compartment within the mouse spleen is comprised of at least five subpopulations, including immature transitional type 1 (T1), transitional type 2 (T2) and T2-preMZ, as well as mature Follicular (Fo) and marginal zone (MZ) B cells (1-4). Signaling through the B cell antigen receptor (BCR) is required for the generation of a functional mature B cell repertoire (3, 5), and multiple components of the BCR signalosome regulate biochemical events after BCR engagement. Upon activation by Btk, PLC- $\gamma$ 2 catalyzes the hydrolysis of phosphatidylinositol bisphosphate (PtdInsP<sub>2</sub> or PIP<sub>2</sub>) into the two second messengers, inositol triphosphate (IP<sub>3</sub>), which activates Ca<sup>2+</sup> signaling, and diacylglycerol (DAG), which activates protein kinase C (PKC) signaling pathways (6, 7). It is well established that both calcium and PKC signaling are important in B cell survival and development/maturation, as well as proliferation (8, 9). Thus, quantities of PIP<sub>2</sub>, DAG and IP<sub>3</sub> are likely to direct the subsequent signaling events leading to the fate of developing B cells. Developing new protocols to study lipid signaling is

therefore a prerequisite for understanding the molecular mechanisms underlying the generation of mature B cells.

Because only a small number of each B cell subpopulation can be obtained by flow cytometry-based cell sorting methods, we developed and optimized mini-scale DAG and IP<sub>3</sub> quantification protocols in order to analyze the levels of DAG and IP<sub>3</sub> within transitional B cell subsets following BCR cross-linking. These protocols allowed the quantification of lipid species utilizing as few as  $0.5-5 \times 10^6$  purified B cells. Additionally, utilizing flow cytometry (FCM) techniques to identify B cell populations, we analyzed PIP<sub>2</sub> levels and calcium mobilization in transitional and mature B cell subsets, again using as few as  $1 \times 10^6$  total B cells. Here we report the detailed experimental procedure of micro-scale flow cytometry and biochemical analysis of lipid signaling in primary B cells.

## MATERIALS AND METHODS

### Mice and B cell purification

Wild type C57BL/6 mice (Jackson Laboratories, Maine, USA) were used. All mice that were used as the source of splenocytes were treated humanely in accordance with federal and state government guidelines and their use was approved by the institutional animal committee.

Primary B lymphocytes were isolated from single cell suspensions pooled from the spleens of C57BL/6 mice. Cell suspensions were depleted of RBCs, and B lymphocytes were purified through magnetic sorting on an AutoMACS automated cell sorter (Miltenyi Biotech) by negative selection (CD43-depletion) to avoid B cell activation. The purity of B cells isolated in this manner was approximately 90–95% as verified by FCM analysis using anti-B220 and anti-IgM antibodies (PharMingen, USA). FACS-sorting for highly purified transitional B cell subsets was performed as previously described (2).

### Lipid analysis


*PIP<sub>2</sub> quantitation:* PIP<sub>2</sub> levels following anti-IgM stimulation were monitored by flow cytometry. MACS-enriched B cells were resuspended ( $10 \times 10^6$  cells/ml) in complete RPMI-1640, rested on ice for 1.5 hours, and incubated for 15 minutes at 37°C prior to assay. Cells were then stimulated with 20 µg/ml goat anti-mouse IgM F(ab')<sub>2</sub> antibodies (Jackson ImmunoResearch) for the indicated times, and immediately fixed (in order to preserve PIP<sub>2</sub> phosphorylation status) by harvesting 100 µl (containing  $1 \times 10^6$  cells) into 1 ml cold PBS + 2% paraformaldehyde at each time-point. Cells were permeabilized utilizing the Cytofix/Cytoperm kit (BD Biosciences) according to the manufacturer's instructions, and intracellularly labeled with a biotinylated antibody raised against PI(4,5)P<sub>2</sub> (Z-B045; Echelon Biosciences, Inc; 1:500 dilution) for 20 minutes on ice. Cells were then washed and labeled with APC-conjugated anti-CD19 (clone 1D3; 1:100 dilution), FITC-conjugated anti-HSA (clone M1/69; 1:500 dilution), and PE-conjugated anti-CD21 (clone 7G6; 1:400 dilution) for the identification of splenic B cell subsets, and PerCP-Cy5.5-conjugated streptavidin (1:300 dilution) to visualize PIP<sub>2</sub>. Anti-C1qRp (AA4.1) and anti-CD23 antibodies were not utilized due to non-reactivity with B cells that have undergone fixation. Following a 20 minute incubation on ice, cells were washed, and resuspended in 100 µl FACS buffer (2% FBS in PBS), followed by 100 µl fix buffer (2% paraformaldehyde in FACS buffer) to preserve the cells. To test the specificity of the anti-PIP<sub>2</sub> antibody, freshly isolated MACS enriched B cells were fixed and permeabilized as described above, then intracellularly labeled with an anti-PIP<sub>2</sub>-FITC antibody (Z-G045; Echelon Biosciences, Inc; 1:300 dilution) that had been pre-incubated with varying concentrations of PIP<sub>2</sub> liposomes (Y-P045 (biotinylated); Echelon Biosciences, Inc) for 15 minutes on ice. Data was acquired on a BD FACSCalibur or

BD LSRII flow cytometer. Results are expressed as MFI and are representative of multiple experiments.

*DAG quantification:* Purified B cells were labeled with [<sup>3</sup>H]palmitic acid (2.5 mCi/ml), as a fatty acid precursor. Labeled B cells were washed twice with PBS, incubated for various time periods or stimulated through the BCR with goat anti-mouse IgM F(ab')<sub>2</sub> as indicated. The reactions were stopped by the addition of excess cold PBS. The cells were pelleted by centrifugation and used for lipid extraction. The pellet was suspended in 2 ml chloroform-methanol (2:1 v/v) and 0.4 ml of 0.5 N HCl was added. The mixture was shaken vigorously, centrifuged, and the upper (aqueous) phase was discarded while the lower (organic) washed twice with theoretical upper phase (10). The radioactivity of lipids was determined on an aliquot fraction by scintillation counting. Dipalmitin carrier standard was then added, the lipid phase evaporated gently under a N<sub>2</sub> stream and the residue dissolved in a minimum volume of chloroform-methanol (2:1 v/v). DAG was separated by thin layer chromatography on silica gel G plates using a double development system and detected under UV light after spraying the plates with primulin as previously described (11). The area containing DAG was scraped and the radioactivity determined as previously described (6, 11, 12).

*[<sup>3</sup>H]-Inositol labeled soluble inositol polyphosphate (IP) profiling:* The soluble IP profile for B cells was determined by adapting a previously described protocol (13). Briefly, MACS-enriched total B cells or FACS-purified B cell subsets ( $10 \times 10^6$ /ml) were cultured in inositol-free DMEM containing 60 mCi [<sup>3</sup>H]Inositol/ml (Perkin Elmer) for 4 hours. Cells were harvested, washed twice with 1 ml PBS and suspended in PBS ( $10 \times 10^6$ /ml). 100 µl cell aliquots (containing  $1 \times 10^6$  cells) were treated with 25 µl PBS with or without 20 µg/ml goat anti-mouse IgM F(ab')<sub>2</sub>. Treatments were stopped by adding 6.6 µl 10 N HCl and chilling cells on ice. Soluble IPs were extracted by adding 465 µl of chloroform/methanol (1:2 v/v). The mixture was vortexed for 2 minutes at maximum speed, followed by addition of 156 µl each of chloroform and 2 M KCl and another 2 minutes of vortexing. The lysate was spun at 13,000 g for 5 minutes and the supernatant recovered. Samples were analyzed by the HPLC (Shimadzu Corporation) connected to the Radio-HPLC Detector (IN/US Systems, Inc.). The IPs were resolved by a Whatman partisphere SAX strong-anion exchange column (4.6 x 125 mm) and a linear gradient from 10 mM to 1.7 M ammonium phosphate (pH 3.5) for 25 minutes, followed by elution with 1.7 M ammonium phosphate for 20 minutes.

*Calcium mobilization in B cell subsets:* MACS-enriched B cells were resuspended ( $10 \times 10^6$  cells/ml) in warm HBSS w/o phenol red, and loaded with indo-1-AM (5 µg/ml) for 30 minutes at 37°C. Cells were then washed, resuspended ( $20 \times 10^6$  cells/ml) in FACS buffer, placed on ice, and labeled for 20 minutes with PE-conjugated anti-HSA (clone M1/69; 1:750 dilution), biotinylated anti-CD23 (clone B3B4; 1:200 dilution) which was visualized by PerCP-Cy5.5-conjugated streptavidin (1:200 dilution), and APC-conjugated anti-C1qRp (clone AA4.1; 1:100 dilution) for B cell subset identification. Antibodies raised against B220,



**Fig. 1A & 1B: PIP<sub>2</sub> is hydrolyzed in all B cell subsets following anti-IgM stimulation.** (A): Freshly isolated MACS-enriched wt B cells were fixed with 2% paraformaldehyde, permeabilized and labeled with antibodies directed against CD19, CD21, and HSA for B cell subset identification. (B): B cells from (A) were examined for basal PIP<sub>2</sub> levels utilizing an anti-PI(4,5)P<sub>2</sub> antibody. Data is an example of multiple experiments.

CD19, and/or CD21 were not used in order to avoid inadvertent stimulation of B cells. Following labeling, cells were washed, resuspended ( $1 \times 10^6$ /ml) in HBSS containing 1 mM calcium, and kept on ice until assay. Prior to acquisition, cells were warmed for 15 minutes in a 37°C water bath, and applied to a BD LSR II cytometer. Baseline, non-stimulated measurements were taken for 90 seconds, followed by the addition of 10  $\mu$ g/ml goat anti-mouse IgM F(ab')<sub>2</sub> antibodies for the indicated time. Ionomycin (1 mM) was used as a positive control to ensure that cells were competent to release calcium from intracellular stores. Data is expressed as a ratio of Indo1-violet/Indo1-blue, and is a representative of multiple experiments.

### Statistics

The data represent the mean of at least three experiments analyzed utilizing Student's *t*-test.

## RESULTS

### PIP<sub>2</sub> hydrolysis occurs in all B cell subsets as determined by FCM

Because PIP<sub>2</sub> serves as a substrate for PLC- $\gamma$ 2 in the production

of DAG (14), we wanted to determine basal and BCR-induced levels of PIP<sub>2</sub> in all B cell subsets. MACS-enriched splenic B cells were stimulated, harvested, and immediately fixed in 2% paraformaldehyde in order to preserve the phosphorylation status of PIP<sub>2</sub>. Cells were then stained intracellular for PIP<sub>2</sub> and extracellular with a combination of anti-HSA, anti-CD21, and anti-CD19 antibodies for B cell subset identification. Utilizing flow cytometric analysis, CD19<sup>+</sup> cells were displayed with HSA and CD21 in order to avoid inclusion of any remaining non-B cells in the analysis (Fig. 1A). CD21<sup>hi</sup> MZ and pre-MZ were separated based on differential HSA expression, which corresponded to CD23 expression (preMZ are HSA<sup>hi</sup>CD23<sup>+</sup>, MZ are HSA<sup>mid</sup>CD23<sup>lo</sup>), as we previously reported (2). Basal PIP<sub>2</sub> levels were then measured in each subset (Fig. 1B). In addition, we determined the kinetics of BCR-induced PIP<sub>2</sub> production in each B cell subpopulation following stimulation with 20  $\mu$ g/ml anti-IgM F(ab')<sub>2</sub> antibodies. In all splenic B cell subsets, PIP<sub>2</sub> levels were reduced immediately in response to stimulation with anti-IgM, indicative of PIP<sub>2</sub> hydrolysis (Fig. 1C).


### Optimal lipid labeling with [<sup>3</sup>H]palmitic acid is achieved in 6 hours

To analyze the lipid synthesis and BCR-directed DAG production in wild type primary B cells, we systematically assessed the numbers of primary B cells and the amount of radioactivity required for optimum labeling. MACS-enriched total splenic B cells or FACS-sorted B cell subsets were used immediately for labeling with [<sup>3</sup>H]palmitic acid, a saturated fatty acid (6). Pilot experiments indicate that optimal lipid synthesis is achieved with a concentration of 2.5 mCi/ml (data not shown). Using this information, we discovered that the incorporation of radioactivity was detectable after 30 min, increased in a time-dependent manner reaching a plateau at 6 h, and remained unchanged through 24 h (Fig. 2A). The same profile was obtained when the results were expressed as cpm/mg protein indicating a dynamic process (data not shown).

### BCR-induced DAG production is dose-dependent and can be detected with as few as $0.5 \times 10^6$ B cells

MACS-enriched total B cells were labeled for 6 h with [<sup>3</sup>H]palmitic acid, then were stimulated for 5 minutes with various concentration of anti-IgM antibodies in order to determine the production of DAG (12). DAG production increases in a dose dependent manner and is optimal at 20  $\mu$ g/ml of anti-IgM (Fab')<sub>2</sub> (Fig. 2B).

To test the feasibility of DAG quantification experiments with small numbers of primary B cells, we labeled MACS-enriched total B cells for 6 h with [<sup>3</sup>H]palmitic acid, and determined its incorporation into whole lipids in various numbers of B cells. We observed a strong correlation between the total radioactivity detected in whole lipids with the number of B cells labeled ( $r = 0.998$ ). It is noteworthy that the radioactivity was close to 50000 and 100000 cpm for  $0.5 \times 10^6$  and  $1 \times 10^6$  labeled cells respectively (Fig. 2C). Since the level of DAG represents 0.6% of total lipids



**Fig. 1C & 1D: PIP<sub>2</sub> is hydrolyzed in all B cell subsets following anti-IgM stimulation. (C):** B cells were stimulated for the indicated times with 20 μg/ml anti-IgM, immediately fixed and permeabilized, then labeled and examined for PIP<sub>2</sub> levels as in (B); T1 (filled diamonds), T2 (filled squares), mature FoB (open circles), preMZ (filled triangles), mature MZB (open triangles). Data is displayed as the average MFI of duplicate samples, and is representative of at least 4 experiments. **(D):** The specificity of the anti-PIP<sub>2</sub> antibody was tested by pre-incubation of the antibody with PIP<sub>2</sub> liposomes for 15 min prior to intracellular labeling of MACS-enriched total B cells as in (B).

extracted from the mouse primary splenic B cells (2), the amount of [<sup>3</sup>H]DAG in 0.5x10<sup>6</sup> B cells should be close to 500 cpm which is more than 100 times greater than the background, and should therefore allow reliable estimation of DAG production.

We next determined the level of BCR-induced DAG production after 5 minutes of stimulation with anti-IgM (20 μg/ml) using various numbers of MACS-enriched total B cells ranging from 0.5x10<sup>6</sup> to 5x10<sup>6</sup>. The results show that upon BCR ligation the profile of DAG production obtained from a small number (0.5x10<sup>6</sup>) of primary B cells coincides precisely with that obtained from a large number (5x10<sup>6</sup>) of B cells (Fig. 2D). When 0.5x10<sup>6</sup> B cells were labeled with [<sup>3</sup>H]palmitic acid, the radioactivity detected in DAG was similar to that calculated from the standard curve (450 cpm corresponds to 0.5% of total lipids). Therefore, lipid signaling in transitional B cell subsets may be investigated using [<sup>3</sup>H]palmitic acid as the radioactive

labeling precursor. In this regard, we previously demonstrated that BCR-mediated DAG production occurs in T2 and mature B cells following both 5 min (2.1 ± 0.1-fold increase, and 5.1 ± 0.3-fold increase, respectively, compared with the non-stimulated control) and 30 min (2.8 ± 0.3-fold increase, and 4.8 ± 0.4-fold increase, respectively, compared with the non-stimulated control) stimulation with anti-IgM (2). In contrast, T1 B cells did not produce DAG in response to anti-IgM stimulation (2).

#### HPLC analysis of BCR-induced IP<sub>3</sub> production in splenic B cell subsets


Because we observed significant differences in BCR-mediated DAG production in transitional B cell subsets, we also investigated the production of IP<sub>3</sub> in response to BCR stimulation utilizing HPLC. Using MACS-enriched total B cells, we first performed a kinetics study and determined that 1 minute of anti-IgM stimulation was the optimum time-point for IP<sub>3</sub> detection (Figure 3A and B). We then utilized 1x10<sup>6</sup> FACS-purified transitional T1 and T2 B cells to investigate IP<sub>3</sub> production following anti-IgM stimulation. T2 B cells produced significantly more IP<sub>3</sub> than T1 B cells (1.52 ± 0.5-fold increase and 3.09 ± 0.2-fold increase, respectively; n = 3) (Figure 3 C-E), consistent with the previously observed DAG production profile (2).

#### FCM analysis of intracellular Ca<sup>2+</sup> mobilization in B cell subsets

Because the generation of IP<sub>3</sub> subsequently induces Ca<sup>2+</sup> mobilization in primary *wt* B cells (15), we next wanted to investigate calcium mobilization in B cell subsets. To precisely define stage-specific calcium mobilization within B cell subsets, we utilized FCM techniques. 1x10<sup>6</sup> MACS-enriched B cells were loaded with indo-1-AM and labeled with a combination of anti-C1qRp (AA4.1), anti-HSA and anti-CD23 antibodies for B cell subset identification (Fig. 4A). Antibodies directed against B220, CD19, and CD21 were not used in order to avoid inadvertent B cell activation. Cells were applied to the cytometer, baseline measurements taken for 60-90 seconds, then stimulated with anti-IgM for the indicated times. From the AA4<sup>+</sup> gate, T1 (HSA<sup>hi</sup>CD23<sup>-</sup>) and T2 (HSA<sup>hi</sup>CD23<sup>+</sup>) cells were gated for calcium analysis. The vast majority of T2-preMZ cells do not express AA4.1; therefore they are excluded from T2 cells gated in this manner. From the AA4<sup>-</sup> gate, mature FoB (HSA<sup>lo</sup>CD23<sup>+</sup>) and MZ (HSA<sup>mid</sup>CD23<sup>lo</sup>) were gated for calcium analysis. Our results demonstrate that calcium mobilization can be monitored in all B cell subsets following anti-IgM stimulation, using as few as 1x10<sup>6</sup> MACS-enriched, FACS-stained total B cells (Fig. 4B).

## DISCUSSION

The goal of this study was to develop protocols to assess analysis of lipid signaling in primary B cell subpopulations




**Fig. 2: Lipid synthesis and DAG production in small numbers of primary B cells.** (A): Kinetics of lipid synthesis in primary B cells. Incorporation of [<sup>3</sup>H]palmitic acid (2.5 mCi/ml) into lipids was determined in purified wt primary B cells. Lipid extraction was performed according the protocol of Folch (10). Lipid synthesis is expressed in cpm and represents the average + SEM of three independent performed in triplicate. (B): BCR-induced DAG production is dose dependent. 5 × 10<sup>6</sup> purified B cells were labeled for 6 h with [<sup>3</sup>H]palmitic acid (2.5 mCi/ml) and stimulated for 5 min with various doses of anti-IgM. Lipids were extracted and separated as described (12) and DAG was quantified. DAG is expressed as percentage increase of DAG of total radioactivity and represents the average + SEM of three independent performed in triplicate. \*P < 0.05; \*\*\*P < 0.001. (C): [<sup>3</sup>H]palmitic acid is incorporated into whole lipids in small numbers of primary B cells. The indicated numbers of non-stimulated B cells from wt mice were labeled for 6 h with [<sup>3</sup>H]palmitic acid (2.5 mCi/ml) and total incorporation of [<sup>3</sup>H]palmitic acid into lipids (r = 0.998) was monitored as described (6). (D): BCR-induced DAG production can be measured in as few as 0.5 × 10<sup>6</sup> primary B cells. The same protocol as in (B) was applied to 0.5–5 × 10<sup>6</sup> primary B cells that were stimulated for 5 minutes with 20 µg/ml of anti-IgM. \*\*\*P < 0.001.

dynamic changes in lipid signaling in each B cell subset following antigenic stimulation. By utilizing FCM techniques to assess both PIP<sub>2</sub> levels and calcium mobilization in individual B cell subsets, we eliminated the necessity for purification of B cell subsets in these analyses. Additionally, we defined optimal conditions to rapidly analyze DAG and IP<sub>3</sub> production utilizing numbers of total primary B cells similar to the quantities of transitional B cells that can be obtained by magnetic bead and FCM based cell sorting protocols (2, 16).

In this report we used a combination of the most current staining schemes to identify B cell subpopulations in the spleen by FCM (1–4, 16, 17). We have combined extracellular staining for B cell subset identification with intracellular staining protocols in order to determine PIP<sub>2</sub> levels in all B cell subsets following stimulation with anti-IgM (Fig. 1). It is important to note that following an initial period of BCR-induced PIP<sub>2</sub> hydrolysis, BCR-induced PIP<sub>2</sub> production occurs in all the B cell subpopulations (Fig. 1C). This increase in PIP<sub>2</sub> over time may result from DAG conversion into PIP<sub>2</sub> or PIP<sub>2</sub> synthesis by PIP5K following BCR stimulation (18, 19). Thus, a PIP<sub>2</sub> cycle is activated in B cells in response to BCR cross-linking. The lower production of DAG and IP<sub>3</sub> in T1 B cells suggests developmental stage-specific mechanisms of PIP<sub>2</sub> cycling in splenic B cell subpopulations. Additionally, it is noteworthy that the B cell subset with the lowest baseline (mature FoB) had the greatest changes in PIP<sub>2</sub> MFI while the B cell subsets with the highest baselines (T2 and preMZ) had the smallest changes in PIP<sub>2</sub> MFI (Fig. 1C), indicating that the inducible production of PIP<sub>2</sub> through BCR stimulation is tightly associated with its basal level. As we have previously shown for PKC phosphorylation (2), this protocol can be extended to study other intracellular biochemical signaling events following antigenic stimulation.

In addition, we have employed FCM techniques in order to assess calcium mobilization in B cell subsets. Again, employing a unique combination of antibodies, which identify B cell subsets but will not activate cells inadvertently, we monitored calcium mobilization in each B cell subset following antigenic stimulation using as few as 1 × 10<sup>6</sup> MACS-enriched B cells and found that all B cell subsets mobilize calcium following anti-IgM stimulation (Fig. 4). These FCM-based techniques use only minimal numbers of primary B cells, and eliminate the need for FACS-based cell sorting of individual B cell populations.

We also report that [<sup>3</sup>H]palmitic acid can be used to efficiently label lipids in primary splenic B cells. We observed that peak radioactivity is reached after 6 h of labeling (Fig. 2A). Since splenic B cells die rapidly *in vitro* in the absence of B cell growth factors and/or agonists, the relatively short duration of labeling (6 h) with [<sup>3</sup>H]palmitic acid minimizes the loss of B cells (16, 17). Our results also show that lipid levels can be determined with small numbers of B cells, and that the radioactivity incorporated into lipids in 0.5 × 10<sup>6</sup> B cells is close to 5 × 10<sup>4</sup> cpm (8000 times higher than the background) which should allow easy quantification of lipid species (Fig. 2C). These data are in agreement with studies performed with primary B cells labeled with [<sup>3</sup>H]arachidonic acid (15), both in




**Fig. 3A & 3B: Generation of IP<sub>3</sub> in primary B cells:** **(A)** Timecourse of IP<sub>3</sub> production in primary B cells.  $3 \times 10^6$  MACS-enriched B cells were labeled with  $^3\text{H}$ inositol (60 mCi/ml) for 4 h, then incubated with PBS (non-stimulated) or stimulated with 20  $\mu\text{g/ml}$  anti-IgM for the indicated times. Production of soluble IP species was measured by HPLC and the position corresponding to each species is indicated. IP<sub>3</sub> production was measured as the total counts detected within the 14.25-15.25 minute time window. HPLC profile of an *ipk1-5* mutant yeast strain (26) serves as reference. **(B):** IP<sub>3</sub> production from the B cells in (A) is expressed as the counts in the position corresponding to IP<sub>3</sub>.

the profile of DAG accumulation and the quantity of radioactivity.

An alternate enzymatic method for the quantification of DAG employed an *Escherichia coli* DAG kinase which catalyzes the production of [ $^{32}\text{P}$ ]phosphatidic acid from DAG substrate (20). This radiosensitive assay was successfully developed for sn-1,2 diacylglycerols quantification. Moreover, the DAG kinase assay was already applied to Epstein-Barr virus negative Burkitt's lymphoma Ramos cells to quantify ceramide (21). Here, our lipid separation method allows the separation of sn-1,2

Antony *et al.* - Micro-scale flow cytometry-based and biochemical analysis of lipid signaling in primary B cell subpopulations  
www.biologicalprocedures.com




**Fig. 3C - 3E: Generation of IP<sub>3</sub> in primary B cells:** **(C):**  $1 \times 10^6$  FACS-purified *wt* T1 and T2 B cells were labeled with  $^3\text{H}$ inositol as in (A) and stimulated with 20  $\mu\text{g/ml}$  anti-IgM for 1 minute. IP<sub>3</sub> production was detected as in (A). Data is a representative of four independent experiments. **(D):** IP<sub>3</sub> production from the T1 and T2 B cells in (C) is expressed as the counts in the position corresponding to IP<sub>3</sub>. **(E)** Fold change in IP<sub>3</sub> level from T1 and T2 B cells following stimulation with anti-IgM was calculated from the counts measured in the position corresponding to IP<sub>3</sub> as in (C and D) ( $n=3$ ).

diacylglycerol from sn-1,3 diacylglycerol species, and the DAG level reflects that of both isomers (12). In addition, we were able to quantify other neutral lipid species (palmitic acid, mTriglycerides, and monoglycerides) in one single experiment (data not shown).

Development of immunocompetent B-lymphocytes is critically dependent on signaling through the BCR (5). Our results show that while 20  $\mu\text{g/ml}$  of anti-IgM (Fab')<sub>2</sub> antibodies induces maximal production of DAG, lower doses also lead to significant production of DAG (Fig. 2B).

We also demonstrate that in response to BCR cross-linking, an increase in DAG production is clearly detectable with as few as  $0.5 \times 10^6$  B cells (Fig. 2D). We have previously employed the methods described in this report to determine DAG production



**Fig. 4: Calcium mobilization in B cell subsets. (A):** MACS-enriched B cells were loaded with Indo-1-AM (5  $\mu\text{g/ml}$ ), and labeled with AA4.1, anti-CD23, and anti-HSA for B cell subset identification. **(B):** Cells were applied to a cytometer, basal  $\text{Ca}^{2+}$  levels were monitored for 90 s, then cells were stimulated with 10  $\mu\text{g/ml}$  anti-IgM (indicated by the arrow) and monitored for 4-6 min; T1 (dotted black line), T2 (black line), mature FoB (gray line), MZ (dotted gray line). Data is a representative of 4-6 independent experiments.

following anti-IgM stimulation in  $1 \times 10^6$  FACS-purified T1, T2 and mature FoB cells (2). In that study, we reported that T2 and FoB cells rapidly produce detectable levels of DAG following anti-IgM stimulation and that increased DAG levels are maintained for a significant period of time, while T1 B cells do not produce DAG at any time-point following anti-IgM stimulation (2).

After we observed that DAG was differentially produced in transitional B cells subsets, we wanted to confirm that  $\text{IP}_3$  production followed a similar pattern, as these two lipids are both by-products of  $\text{PIP}_2$  hydrolysis. Consistent with DAG production,  $\text{IP}_3$  was detected by HPLC in total B cells and in purified T2 B cells, but not in purified T1 B cells (Fig. 3 A and C). In addition to  $\text{IP}_3$ , we found that levels of other IPs (i.e.,  $\text{IP}_1$ ,  $\text{IP}_2$ ,  $\text{IP}_4$ , and  $\text{IP}_5$ ) also change in all B cell populations, however, these changes are most apparent in total B cells (Fig. 3A). This observation may reflect more robust changes in mature follicular B cells, which is the predominant population Antony *et al.* - Micro-scale flow cytometry-based and biochemical analysis of lipid signaling in primary B cell subpopulations [www.biologicalprocedures.com](http://www.biologicalprocedures.com)

(approximately 60-70%) in splenic B cells. Thus, the magnitude of BCR-induced  $\text{IP}_3$  production is likely to be underestimated because  $\text{IP}_3$  can undergo additional phosphorylations by kinases to generate  $\text{IP}_4$ ,  $\text{IP}_5$ ,  $\text{IP}_6$ . These findings also raise the possibility that other IPs may also play distinct roles in modulating B cell physiology. Indeed, a previous report showed that the PH domain of Btk functions as an inositol 1,3,4,5-tetrakisphosphate ( $\text{IP}_4$ ), inositol 1,3,4,5,6-pentakisphosphate ( $\text{IP}_5$ ), and inositol 1,2,3,4,5,6-hexakisphosphate ( $\text{IP}_6$ ) binding domain and that *xid* mutations (R28C) found in the Pleckstrin Homology domain result in a dramatic reduction of  $\text{IP}_4$  binding activity (22). Consistent with these results, BCR-induced DAG/ $\text{IP}_3$  production is Btk dependent, suggesting that B cell differentiation is likely to be correlated with the  $\text{IP}_4$  binding capacity of the PH domain of Btk (2, 22, 23). Additionally, a recent report revealed that inositol 1,4,5-triphosphate 3-kinase and its reaction product  $\text{IP}_4$  are both required to mediate survival signalling in B cells via a Ras-ERK signalling pathway that controls the expression of the proapoptotic protein, Bim (24).

Together, these results demonstrate our ability to perform experiments to determine BCR-dependent lipid signaling utilizing the small numbers of B cells that are obtained with any of the B cell subpopulation purification protocols we have previously published (2, 16). Additionally, while these protocols were established in cells that were stimulated through the BCR, it is reasonable to believe that these protocols will be useful when monitoring B cell activation by many different cellular stimuli; in fact, we have already utilized these protocols to monitor DAG and  $\text{IP}_3$  production, as well as calcium mobilization, in B cells stimulated with BAFF (25). Thus, the protocols that we describe in this report should allow a comprehensive determination of differential lipid signaling and signal transduction in primary B cell populations in response to various stimuli.

## ACKNOWLEDGMENTS

This study was supported in part by the NIH grant R01 AI060729 (to W.N.K.). The authors wish to thank John Lowe and Joan M. Llanes for technical assistance, and Jim Higginbotham and Michele Nadeff of the VUMC Flow Cytometry Core for assistance with Calcium experiments.

## REFERENCES

- Allman D, Lindsley RC, DeMuth W, Rudd K, Shinton SA, Hardy RR. Resolution of three nonproliferative immature splenic B cell subsets reveals multiple selection points during peripheral B cell maturation. *J Immunol* 2001;167(12):6834-40.
- Hoek KL, Antony P, Lowe J, Shinnars N, Sarmah B, Wentz SR, et al. Transitional B cell fate is associated with developmental stage-specific regulation of diacylglycerol and calcium signaling upon B cell receptor engagement. *J Immunol* 2006;177(8):5405-13.
- Loder F, Mutschler B, Ray RJ, Paige CJ, Sideras P, Torres R,

- et al. B cell development in the spleen takes place in discrete steps and is determined by the quality of B cell receptor-derived signals. *J Exp Med* 1999;190(1):75-89.
4. Srivastava B, Quinn WJ, 3rd, Hazard K, Erikson J, Allman D. Characterization of marginal zone B cell precursors. *J Exp Med* 2005;202(9):1225-34.
  5. Su TT, Guo B, Wei B, Braun J, Rawlings DJ. Signaling in transitional type 2 B cells is critical for peripheral B-cell development. *Immunol Rev* 2004;197:161-78.
  6. Antony P, Petro JB, Carlesso G, Shinnars NP, Lowe J, Khan WN. B cell receptor directs the activation of NFAT and NF-kappaB via distinct molecular mechanisms. *Exp Cell Res* 2003;291(1):11-24.
  7. Coggeshall KM, Cambier JC. B cell activation. VIII. Membrane immunoglobulins transduce signals via activation of phosphatidylinositol hydrolysis. *J Immunol* 1984;133(6):3382-6.
  8. Hikida M, Johmura S, Hashimoto A, Takezaki M, Kurosaki T. Coupling between B cell receptor and phospholipase C-gamma2 is essential for mature B cell development. *J Exp Med* 2003;198(4):581-9.
  9. Wen R, Chen Y, Xue L, Schuman J, Yang S, Morris SW, et al. Phospholipase Cgamma2 provides survival signals via Bcl2 and A1 in different subpopulations of B cells. *J Biol Chem* 2003;278(44):43654-62.
  10. Folch J, Lees M, Sloane Stanley GH. A simple method for the isolation and purification of total lipides from animal tissues. *J Biol Chem* 1957;226(1):497-509.
  11. Antony P, Farooqui AA, Horrocks LA, Freysz L. Effect of D609 on phosphatidylcholine metabolism in the nuclei of LA-N-1 neuroblastoma cells: a key role for diacylglycerol. *FEBS Lett* 2001;509(1):115-8.
  12. Antony P, Petro JB, Carlesso G, Shinnars NP, Lowe J, Khan WN. B-cell antigen receptor activates transcription factors NFAT (nuclear factor of activated T-cells) and NF-kappaB (nuclear factor kappaB) via a mechanism that involves diacylglycerol. *Biochem Soc Trans* 2004;32(Pt 1):113-5.
  13. York JD, Odom AR, Murphy R, Ives EB, Wentz SR. A phospholipase C-dependent inositol polyphosphate kinase pathway required for efficient messenger RNA export. *Science* 1999;285(5424):96-100.
  14. Chien MM, Cambier JC. Divalent cation regulation of phosphoinositide metabolism. Naturally occurring B lymphoblasts contain a Mg2(+)-regulated phosphatidylinositol-specific phospholipase C. *J Biol Chem* 1990;265(16):9201-7.
  15. Bijsterbosch MK, Meade CJ, Turner GA, Klaus GG. B lymphocyte receptors and polyphosphoinositide degradation. *Cell* 1985;41(3):999-1006.
  16. Petro JB, Gerstein RM, Lowe J, Carter RS, Shinnars N, Khan WN. Transitional type 1 and 2 B lymphocyte subsets are differentially responsive to antigen receptor signaling. *J Biol Chem* 2002;277(50):48009-19.
  17. Su TT, Rawlings DJ. Transitional B lymphocyte subsets operate as distinct checkpoints in murine splenic B cell development. *J Immunol* 2002;168(5):2101-10.
  18. Saito K, Tolias KF, Saci A, Koon HB, Humphries LA, Scharenberg A, et al. BTK regulates PtdIns-4,5-P2 synthesis: importance for calcium signaling and PI3K activity. *Immunity* 2003;19(5):669-78.
  19. Toker A. The synthesis and cellular roles of phosphatidylinositol 4,5-bisphosphate. *Curr Opin Cell Biol* 1998;10(2):254-61.
  20. Preiss J, Loomis CR, Bishop WR, Stein R, Niedel JE, Bell RM. Quantitative measurement of sn-1,2-diacylglycerols present in platelets, hepatocytes, and ras- and sis-transformed normal rat kidney cells. *J Biol Chem* 1986;261(19):8597-600.
  21. Kroesen BJ, Pettus B, Luberto C, Busman M, Sietsma H, de Leij L, et al. Induction of apoptosis through B-cell receptor cross-linking occurs via de novo generated C16-ceramide and involves mitochondria. *J Biol Chem* 2001;276(17):13606-14.
  22. Fukuda M, Kojima T, Kabayama H, Mikoshiba K. Mutation of the pleckstrin homology domain of Bruton's tyrosine kinase in immunodeficiency impaired inositol 1,3,4,5-tetrakisphosphate binding capacity. *J Biol Chem* 1996;271(48):30303-6.
  23. Khan WN, Alt FW, Gerstein RM, Malynn BA, Larsson I, Rathbun G, et al. Defective B cell development and function in Btk-deficient mice. *Immunity* 1995;3(3):283-99.
  24. Marechal Y, Pesesse X, Jia Y, Pouillon V, Perez-Morga D, Daniel J, et al. Inositol 1,3,4,5-tetrakisphosphate controls proapoptotic Bim gene expression and survival in B cells. *Proc Natl Acad Sci U S A* 2007;104(35):13978-13983.
  25. Shinnars NP, Carlesso G, Castro I, Hoek KL, Corn RA, Woodland RL, et al. Bruton's Tyrosine Kinase Mediates NF-kappaB Activation and B Cell Survival by B Cell-Activating Factor Receptor of the TNF-R Family. *J Immunol* 2007;179(6):3872-80.
  26. Ives EB, Nichols J, Wentz SR, York JD. Biochemical and functional characterization of inositol 1,3,4,5, 6-pentakisphosphate 2-kinases. *J Biol Chem* 2000;275(47):36575-83.
- Antony *et al.* - Micro-scale flow cytometry-based and biochemical analysis of lipid signaling in primary B cell subpopulations  
www.biologicalprocedures.com

## PROTOCOLS

### B cell purification by negative selection

#### Materials:

FACS buffer (PBS + 2% FBS)  
 MACS buffer (PBS + 0.5% BSA and 1mM EDTA)  
 Complete RPMI (RPMI + 10% FBS, 4mM L-glut, 100U P/S, 0.1% 2-ME, 10mM HEPES)  
 \* anti-mouse CD43 magnetic beads (130-049-801 Miltenyi Biotech)  
 \* AutoMACS automated cell sorter or MiniMACS columns (Miltenyi Biotech)  
 \*\* B-cell enrichment kit (557792 BD Biosciences) and manual IMAG magnet

- 1) Harvest the spleen from a wt C57BL/6 mouse.
- 2) Isolate splenocytes by smashing spleen in a 70 mM cell strainer with 5ml PBS. Transfer cells to a 15ml conical tube and rinse strainer 2X with 5ml PBS. Transfer washes to the tube with cells and centrifuge for 6 minutes at 1600 rpm at 4°C.
- 3) Aspirate supernatant and lyse RBC, while vortexing, by the addition of 900  $\mu$ L sterile filtered ddH<sub>2</sub>O followed immediately by the addition of 100  $\mu$ L 10x PBS. Remove floating fat, and add FACS buffer to 10 ml.
- 4) Count cells.
- 5) Centrifuge cells for 6 minutes at 1600 rpm at 4°C. Resuspend cells at 100x10<sup>6</sup>/ml in MACS buffer.
- 6) \*Add anti-CD43 beads and incubate according to the manufacturers instructions.
- 7) \*Sort cells on the AutoMACS using the DelpeteS cycle or on the MiniMACS magnetic columns according to the manufacturers directions.
- 8) Resuspend CD43-negative cells in 1-2 ml of complete RPMI and count cells.
- 9) Stain a small aliquot for purity with anti-B220 or anti-CD19, and analyze by FCM.

\*\* 6-7) Alternatively, B cells can be enriched using the B cell enrichment kit and manual magnet from BD according to the manufacturers instructions.

### PIP<sub>2</sub> quantitation by FCM

#### Materials:

Stop buffer (PBS + 2% paraformaldehyde)  
 BD Cytotfix/Cytoperm kit  
 Fix buffer (FACS + 2% paraformaldehyde)  
 \*Biotinylated anti-PIP<sub>2</sub> antibody (Z-B045; Echelon Biosciences)  
 Anti-HSA-FITC antibody (BD clone M1/69)  
 Anti-CD21-PE antibody (BD clone 7G6)  
 \*SAv-PerCP-Cy5.5 (BD 551419)  
 Anti-CD19-APC antibody (BD clone 1D3)  
 \*anti-PIP<sub>2</sub>-FITC (Z-G045; Echelon Biosciences)

- 1) Resuspend MACS-enriched B cells in RPMI (10x10<sup>6</sup>/ml) and rest on ice for 1hr.
- 2) Warm cells in a 37°C water bath for 15 minutes
- 3) Prepare "Stop" tubes (1ml stop buffer) for each time point
- 4) Harvest Time 0 sample by collecting 100  $\mu$ l (1x10<sup>6</sup> cells) into the appropriate stop tube (harvest an extra sample if staining with an isotype control)
- 5) Stimulate cells for desired time with anti-IgM F(ab')<sub>2</sub>
- 6) Harvest cells by collecting 100  $\mu$ l at each desired time point into appropriate stop tubes
- 7) Centrifuge 1700 rpm for 4 minutes
- 8) Fix/permeabilize using BD Cytotfix/cytoperm kit according to manufacturer's instructions
- 9) \*Intracellular stain with anti-PIP<sub>2</sub>-biotin (1:500) in 100  $\mu$ l perm-wash buffer/1x10<sup>6</sup> cells for 20 minutes. Wash with 1ml perm-wash, and spin at 1700 rpm for 4 minutes.
- 10) \*Stain intracellular with SAv-PerCP-Cy5.5 (1:300), and extracellular with anti-HSA (1:500), anti-CD21 (1:250), and anti-CD19 (1:100) in 100  $\mu$ l perm-wash buffer/1x10<sup>6</sup> cells for 20 minutes. Wash with 1ml perm-wash, and spin at 1700 rpm for 4 minutes.
- 11) Resuspend cells in 100 $\mu$ L FACS + 100 ml Fix.

12) Acquire data on FACS machine.

\*\* 9-10) Alternatively, anti-PIP<sub>2</sub>-FITC (1:300) can be used to study PIP<sub>2</sub> levels in total B cells (for subset identification, use in combination with HSA-PE (1:500), CD19-PerCP-Cy5.5 (1:200), and CD21-APC (1:400); all from BD)

### DAG quantification by TLC

- 1) Resuspend purified B cells in complete RPMI, and label cells with 2.5 mCi/ml [<sup>3</sup>H]palmitic acid (Sigma) for 6 hours at 37°C, 5% CO<sub>2</sub>.
- 2) Spin cells at 2000 rpm for 5 minutes.
- 3) Remove supernatant
- 4) Wash cells twice with RPMI without serum
- 5) Resuspend cells (3x10<sup>6</sup>/ml) in PBS in glass tubes and warm to 37°C .
- 6) Incubate cells with 20 µg/ml anti-IgM F(ab')<sub>2</sub> for appropriate times at 37°C.
- 7) Stop reactions by adding 1 ml cold PBS
- 8) Spin the cells at 3000 rpm for 5 minutes and wash twice with PBS
- 9) Remove supernatant
- 10) Extract lipids with 2 ml of a mixture of cold chloroform:methanol (2:1, v:v)
- 11) Place the samples at 4 °C overnight
- 12) Add 400 µl HCl (0.5N) to the 2 ml sample
- 13) Vortex tubes for 10 seconds
- 14) Spin at 3000 rpm for 5 minutes at 4°C
- 15) Take 10 µl of the upper phase and count the radioactivity in scintillation counter to ensure that after washing the cells no radioactivity remains in the aqueous phase. Presence of radioactivity in the aqueous phase will indicate poor lipid labelling and the experiment should be discarded.
- 16) Remove the upper phase
- 17) Transfer 1 ml of the organic phase to a new glass tube
- 18) Remove 10 µl of the organic phase and count the radioactivity in scintillation counter
- 19) Add 40 µl of cold carrier lipids (mixture of phosphatidic acid, monoglycerides, 1-2 dipalmitin, palmitic acid, tripalmitin, each at 10 mg/ml; Sigma)
- 20) Evaporate the organic phase gently under an N<sub>2</sub> stream
- 21) Resuspend the residue in 40 µl chloroform-methanol (2:1, v:v)
- 22) Load on a 20x20 cm Silica Gel plate (Merck)
- 23) Separate neutral lipids from the other lipids in the upper phase with a mixture of ethylacetate:isooctane:acetic acid:water (130:20:30:100)
- 24) Allow the TLC to run half-way on the plate (25-30 minutes) in the first dimension
- 25) Dry the plate gently with a hair dryer for 10 minutes
- 26) Further separate the neutral lipids with a mixture of isopropylether:acetic acid (96:4, v:v) in the second dimension
- 27) Dry the plates gently with a hair dryer for 10 minutes
- 28) Spray the plates with primulin (in acetone-water 8:2, v:v)
- 29) Using the carrier lipids as markers, identify the lipid species, and scrape the desired lipids (in silica) into scintillation vials
- 30) Add 150 µl of water and 3 ml of scintillation fluid
- 31) Vortex each tube for 20 seconds
- 32) Count radioactivity

### IP<sub>3</sub> quantification by HPLC

#### *Soluble Inositol Polyphosphate (IP) Extraction*

- 1) Incubate MACS-enriched total B cells or FACS-purified B cell subsets (10x10<sup>6</sup>/ml) in inositol-free DMEM containing 60 mCi [<sup>3</sup>H]inositol/ml (Perkin Elmer) for 4 hours.
- 2) Harvest cells, wash twice with 1 ml PBS and suspended in PBS (10x10<sup>6</sup>/ml) in a glass tube.
- 3) Dispense 100 µl aliquots (1x10<sup>6</sup> cells) in 2.0 ml screw-cap polypropylene tubes.
- 4) Warm cells to 37°C for 5 min. Treat cells with 25 µl PBS with or without 20 µg/ml goat anti-mouse IgM F(ab')<sub>2</sub> for desired time.
- 5) Stop treatments by adding 6.6 µl 10 N HCl and chilling cells instantly on ice.

- 6) Extract soluble IPs by adding 465  $\mu$ l of chloroform/methanol (1:2 v/v). Vortex the mixture at maximum speed for 2 minutes.
- 7) Add 156  $\mu$ l each of chloroform and 2 M KCl and vortex for another 2 minutes.
- 8) Spin the lysate at 13,000 g for 5 minutes and collect the supernatant in a screw-cap tube. Store at -20°C if needed.

### Separation of IPs through HPLC

Analyze the samples by the HPLC (Shimadzu Corporation) connected to the Radio-HPLC Detector (IN/US Systems, Inc.).

Sample preparation: Add ammonium phosphate (pH 3.5) to the supernatant (10 mM final). Filter through an Amicon Ultrafree-MC 0.45 mm filter unit (Millipore).

HPLC: Inject sample (500  $\mu$ l) onto a Whatman partisphere SAX strong-anion exchange column (4.6 x 125 mm). Resolve the IPs by a linear gradient from 10 mM to 1.7 M ammonium phosphate (pH 3.5) for 25 minutes, followed by elution with 1.7 M ammonium phosphate for 20 minutes. For detection by an online Radio-HPLC detector, use Tru-Count (IN/US System) scintillation counting solution in 1:1 ratio.

### Calcium mobilization by FCM

Materials needed:

- HBSS w/o phenol red (containing 10% FCS, 1mM Mg<sup>++</sup>, with/without 1mM Ca<sup>++</sup>)
- Indo-1-AM (Molecular probes I-1223), resuspended at 1 mg/ml in DMSO
- Anti-HSA-PE antibody (BD clone M1/69)
- Anti-CD23-biotinylated antibody (BD clone B3B4)
- AA4.1-APC (e-bioscience 17-5892)
- Non-indo-labeled cells to use for single stain controls for FCM

- 1) Rest MACS-enriched B cells on ice in complete RPMI for 1-2hrs
- 2) Wash cells, resuspend ( $10 \times 10^6$  cells/ml) in warm HBSS (without calcium), and load with indo-1-AM (5 $\mu$ g/ml) for 30 minutes in a humidified 37°C incubator.
- 3) Wash cells twice with HBSS
- 4) Suspend cells ( $20 \times 10^6$  cells/ml) in FACS, and place on ice (reserve a small aliquot of indo-labeled cells to serve as the indo-1 single stain control for FCM)
- 5) Label remaining cells with biotinylated-CD23 (1:200 dilution) for 20 minutes (label single stain control cells as well)
- 6) Wash cells, and suspend ( $20 \times 10^6$  cells/ml) in FACS
- 7) Label cells with HSA-PE (1:750 dilution), SA<sub>v</sub>-PerCP-Cy5.5 (1:200 dilution), and AA4-APC (1:100 dilution) for 20 minutes (label single stain control cells as well)
- 8) Wash cells, and suspend ( $1 \times 10^6$  cells/ml) in HBSS (with calcium)
- 9) Keep on ice until assay
- 10) Warm cells for 15 minutes in a 37°C water bath prior to assay
- 11) Compensate the flow cytometer with single stain control cells
- 12) Load sample onto flow cytometer set to measure the ratio of Indo-1-violet/Indo-1-blue
- 13) Take baseline (non-stimulated) measurement for 30-90 seconds
- 14) Stimulate cells, and measure calcium ratio for 4-10 minutes

\*\* To monitor calcium efflux from the ER without influx of extracellular calcium, use calcium-free HBSS when acquiring. Calcium can be added back in to monitor calcium influx.